

Fiji workshop 2025

- Day 1:
 - Image data visualization in Fiji
- Day 2:
 - Quantitative image analysis in Fiji
- Day 3:
 - One-on-one analysis support with your own image data

Joel Ryan, Philip Kesner, Barbara Da Rocha, Mar Garcia Ferrés

ABIF


March 4 to 6 2025

<https://github.com/ABIF-McGill/ABIF-Fiji-Workshop-2025>

Day 1: Data visualization in Fiji

- Today's goals
 - Get acquainted with Fiji!
 - LUTs, contrast
 - Composite images
 - Intro to scripting!
- Look at data from 3D images
- Look at data from timelapse images
- Images for publications


Fiji is awesome! But it has some drawbacks...



- Why should you use Fiji?
 - Free open-source software
 - Huge development community
 - Libraries to open nearly any microscope image (!)
 - Automation and scripting with multiple programming languages (!)
- When should you use something else than Fiji?
 - More powerful computation
 - Better 3D visualization
 - Many deep learning tools
 - Instrument-specific processing

Digital images

- A digital image is simply a table of numbers
- Each number represents the intensity value of a pixel
- The intensity value is proportional to the amount of light collected at that location
- Software such as Fiji simply displays each number as a certain shade or colour



Example images






Figure 1: Images of CHO-K1 cells stained with Phalloidin Alexa Fluorophore 488 and stained with the nuclear probe DAPI. Images were flat field corrected using images of a fluorescent plastic slide imaged on the same microscope. Images were collected on an upright Zeiss Axioskop with a EC PlanNeoFluar 20x/0.5 NA objective lens with an AxioCam ICm1 camera. DAPI was imaged using a DAPI cube and an exposure time of 50 ms. Actin was imaged with a FITC cube and a 150 ms exposure. Brightness, gamma and contrast were adjusted to visualize the actin features well. The scale bar is 50 µm. Box plot of actin intensity for control (n=110) and treated (n=149) cells. T-test was done with a two sample t-test, unequal variance, two tailed. *** corresponds to $P < 0.001$.


Example images




C. elegans embryo
eGFP-H2B
eGFP-gamma-tubulin
Max projection movie
One z-stack per minute acquisition
4 frames-per-second playback




0 minutes



2 minutes




4 minutes




8 minutes

Example images



Imaris demo image

Pyramidal cell z-stack
69 z-slices
15 frames-per-second playback






Maximum intensity projection (Fiji)
Scalebar: 10 um

Example images




Imaris demo image

Drosophila egg Chamber
3 channels x 98 z-slices
16 frames-per-second playback




Orthogonal slice views

Exercise 1 – raw vs composite RGB images



demo_DAPI_Phalloidin_Mitotracker_001.tif

Raw image, drag-and-drop in powerpoint.



demo_DAPI_Phalloidin_Mitotracker_001.tif (RGB).tif


Composite RGB image generated with Exercise 1 macro,
drag-and-drop in powerpoint.

Day 2: Quantitative imaging analysis in Fiji


- Goal:
 - Learn to use common image processing tools in Fiji in order to extract quantitative information

Introduction to quantitative imaging analysis


- Why quantitative analysis?



Tet1




Tet1 $\Delta 1-131$




Tet2

Introduction to quantitative imaging analysis


- Why quantitative analysis?




Tet1



Tet1 $\Delta 1-131$




Tet2




Introduction to quantitative imaging analysis


- Why quantitative analysis?




Tet1



Tet1 Δ 1-131



Tet2



- Measured quantity (space, time, intensity etc)
- Increases robustness of observations
- Elucidate patterns that are not easily visible by eye




Better questions, better science!

Introduction to quantitative imaging analysis

Images can be quantitative:

- A digital image is simply a table of numbers
- Each number represents the intensity value of a pixel
- The intensity value is proportional to the amount of light collected at that location




Value at each pixel = **Signal + Background + Noise**

Introduction to quantitative imaging analysis

- Goal:
 - Learn to use common image processing tools in Fiji in order to extract quantitative information

Introduction to quantitative imaging analysis


- Goal:
 - Learn to use common image processing tools in Fiji in order to extract quantitative information




Open an image

Introduction to quantitative imaging analysis

- Goal:
 - Learn to use common image processing tools in Fiji in order to extract quantitative information




Open an image




Processing such
as filtering

Introduction to quantitative imaging analysis


- Goal:
 - Learn to use common image processing tools in Fiji in order to extract quantitative information



Open an image




Processing such
as filtering




Segmentation or
object detection

Introduction to quantitative imaging analysis

- Goal:
 - Learn to use common image processing tools in Fiji in order to extract quantitative information




Open an image



Processing such as filtering



Segmentation or object detection



Data extraction from segmented objects

- Number of objects
- Shape / size
- Intensity within object...

Introduction to quantitative imaging analysis

- Goal:
 - Learn to use common image processing tools in **Fiji** in order to extract quantitative information




- Why Fiji?
 - Free open-source software
 - Huge development community
 - Libraries to open nearly any microscope image (!)
 - Automation and scripting with multiple programming languages (!)

Introduction to quantitative imaging analysis

- **Today:**

- Load images in Fiji, get familiarised with the program
- Processing using convolutional filters
- Segmentation with thresholding
- Data extraction
 - ...plotting?



for loop

```
for (i = 0; i < numROI; i++) {  
  
    roiManager("Select", i);  
    Stack.setChannel(2);  
    run("Measure");  
  
}
```


for loop

Start value

End value

Increment

```
for (i = 0; i < numROI; i++) {  
  
    roiManager("Select", i);  
    Stack.setChannel(2);  
    run("Measure");  
  
}
```

for loop

Start value

End value

Increment

```
for (i = 0; i < numROI; i++) {  
  
    roiManager("Select", i);  
    Stack.setChannel(2);  
    run("Measure");  
  
}
```

first, “build a sequence of integers going from the Start Value, to the End value, with an increment of 1”

for loop

Start value

End value

Increment

```
for (i = 0; i < numROI; i++) {  
  
    roiManager("Select", i);  
    Stack.setChannel(2);  
    run("Measure");  
  
}
```

first, “build a sequence of integers going from the Start Value, to the End value, with an increment of 1”

---in the case of numROI being 5:

0, 1, 2, 3, 4

for loop

Start value

End value

Increment

```
for (i = 0; i < numROI; i++) {  
  
    roiManager("Select", i);  
    Stack.setChannel(2);  
    run("Measure");  
  
}
```

first, “build a sequence of integers going from the Start Value, to the End value, with an increment of 1”

---in the case of numROI being 5:

0, 1, 2, 3, 4,

“Run the { content of the loop } replacing i with the first integer of the sequence.

Then run the { content of the loop } replacing i with the second integer in the sequence... “ and so on...

for loop

Start value

End value

Increment

```
for (i = 0; i < numROI; i++) {
```

```
    roiManager("Select", i);
```

```
    Stack.setChannel(2);
```

```
    run("Measure");
```

```
}
```

```
i = 0  
roiManager("Select", i);  
Stack.setChannel(2);  
run("Measure");
```

```
i = 1  
roiManager("Select", i);  
Stack.setChannel(2);  
run("Measure");
```


```
i = 2  
roiManager("Select", i);  
Stack.setChannel(2);  
run("Measure");
```

```
i = 3  
roiManager("Select", i);  
Stack.setChannel(2);  
run("Measure");
```


```
...
```


Deep Learning based segmentation

- Sometimes objects are too close for adequate thresholding




DAPI image



Otsu threshold image

Deep Learning based segmentation


- Machine learning tools trained on cell images can help distinguish close objects



DAPI image




Nuclei masks identified with
Cellpose




DAPI image + mask outlines

Deep Learning based segmentation


- Using the DAPI masks, we can then quantify intensities of other stainings in the nucleus



DAPI image + nucleus mask outlines




LIN28-GFP image



LIN28-GFP image + nucleus mask outlines


Deep Learning based segmentation

- Using the DAPI masks, we can then quantify intensities of other stainings in the nucleus



LIN28-GFP image + nucleus mask outlines

Results					Area	Mean	StdDev	Mode
Label								
342	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0949-0908			147.453	12.157	5.874	10
343	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0944-0376			155.523	32.527	17.836	20
344	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0935-0813			66.584	12.790	7.279	10
345	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0950-0581			187.322	26.753	14.969	20
346	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0951-0448			148.421	36.044	22.070	20
347	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0952-1002			116.864	8.281	4.259	5
348	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0960-0678			87.810	22.566	14.736	11
349	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0963-0029			133.087	27.723	14.907	15
350	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0963-0517			83.694	29.380	21.045	12
351	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0971-0955			132.118	7.816	3.847	5
352	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0972-0146			127.518	30.303	17.061	20
353	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0972-0813			83.774	16.490	11.709	7
354	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0959-0852			25.905	13.484	0.645	8
355	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0967-0076			44.793	60.728	21.658	66
356	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0991-0640			192.165	28.741	13.635	19
357	MAX_20240129_CGR8_iBidi_Timecourse_OTX2_SOX2_ESRRB_LIN28A_CGR8_OTX2_T0_1-2.tif	0994-0918			115.089	12.575	5.527	8



SuperPlots: Communicating reproducibility and variability in cell biology

- SuperPlots provides guidelines for presenting experimental data.
- Best to show all data points and colour code them per experimental replicate.
- For statistics the n value should be the number of experiments.




Plotting resources

PlotsOfData - Plots all Of the Data

Data upload

- ☒ Example 1 (wide format)
- ☐ Example 2 (tidy format)
- ☐ Upload file
- ☐ Paste data
- ☐ URL (csv files only)
- ☐ These data are Tidy

Select and order:

 Download in tidy format (csv)

☐ Show information on data formats

Data offset

- ☒ Quasirandom
- ☐ Random
- ☐ None; stripes
- ☐ None (for small n)

Visibility of the data


0

0.3

1

Statistics

- ☒ Median
- ☐ Mean
- ☐ Boxplot (minimum n=10)
- ☐ Violin Plot (minimum n=10)
- ☐ Add 95% CI (minimum n=10)



Generate plots with data points on PlotsOfData webapp

<https://huygens.science.uva.nl/PlotsOfData/>